ULTRAFILTER

U 8000 S

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 **GAMBRO**

General Information

The information contained in this instruction manual must be read before the use of this device to ensure a safe and effective treatment. Review all WARNINGS and CAUTIONS before use and refer to specific operational WARNINGS and CAUTIONS described in the text whenever necessary.

The U 8000 S is delivered sterile with non-pyrogenic fluid compartments.

The quality of the product is guaranteed only when it is prepared and used according to the procedures described in the "Directions for use".

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

WARNING ! is used to alert the operator not to take a specific action that can cause potential hazard which, if ignored, can result in an adverse reaction, injury or death.

CAUTION ! is used to alert the operator to take a specific action to protect against possible hazard which, if ignored, could have an adverse effect on the patient or equipment.

NOTE ! is used to remind the operator of normal treatment functions and what is the suitable action in a particular situation.

INDICATIONS

Ultrafiltration with a Gambro ultrafilter is indicated for purification of incoming water for dialysis as well as purification of dialysis fluid to obtain microbiologically high quality dialysis fluid. The risk of exposure to bacteria and endotoxins can thereby be minimized. In combination with an additional single use ultrafilter the U 8000 S can be also used for the on-line preparation of non-pyrogenic substitution fluid with the Gambro AK 100/200 ULTRA and AK 200 ULTRA S.

CONTRAINDICATIONS

None known.

ADVERSE REACTIONS

If the dialysis fluid or the incoming water is heavily contaminated or the handling of the U 8000 S is incorrect, there is a risk of bacterial and/or pyrogenic overload in the ultrafilter with consequent risk of breakthrough. The patient may then be exposed to bacteria or pyrogens sufficient to induce a pyrogenic reaction with symptoms of shivering and raised body temperature. Pyrogenic reaction may be accompanied by severe reduction in blood pressure.

Patient exposure to high concentrations of disinfectants due to inadequate rinsing after chemical disinfection of the U 8000 S may lead to symptoms of rash, chest pain and hemolysis. Exposure to formaldehyde may lead to the development of anti-N-like antibodies with risk of hemolytic anemia. It is, therefore, imperative to validate the rinsing procedure after each disinfection procedure (Residual test).

In the event of a severe patient reaction, the blood lines shall be clamped immediately and appropriate medical treatment initiated. The blood in the hemodialyzer and blood lines should not be returned to the patient and it is recommended that the treatment session be discontinued. The U 8000 S shall be discarded and the dialysis machine disinfected. A new sterile U 8000 S can then be inserted in the dialysis fluid circuit according to instructions under DIRECTIONS FOR USE.

GENERAL WARNINGS AND CAUTIONS

The U 8000 S has been designed for the use with Gambro dialysis machines. If other dialysis machines are used it has to be ensured that they are able to compensate for the increased pressure drop caused by the insertion of the U 8000 S in the fluid circuit.

The user must ensure that the increased pressure drop does not interfere with the control of ultrafiltration (the patient's weight loss) during treatment. Note that the increased pressure drop may also limit the use of certain types of dialyzers due to specific pressure regulations in the machine.

Various functions of the dialysis machine may also be disturbed causing alarm situations.

Make certain that the procedures outlined in this package insert are carefully studied before the U 8000 S is used.

An aseptic technique must be employed to avoid contamination of the fluid pathways.

The U 8000 S is designed for multiple use. Gambro can only guarantee the performance, sterility and non-pyrogenicity of the U 8000 S if it has been sterilized by Gambro or by a subcontractor approved by Gambro.

It is the sole responsibility of the physician in charge to validate the microbiological quality of the system in order to prove that the applied disinfection process is effective and that the retention performance of the filter is still appropriate. A regular monitoring of both water and filtrate quality is advised. The ultrafilters used as water filter and dialysis fluid filter are to be changed routinely after use for one to two months maximum, depending on the result of the microbiological controls.

It is imperative to use only chemical disinfectants recommended in this operator's manual; see section "Directions for use": Treatment. The instructions supplied by the manufacturer of the chemical disinfectant agent must be strictly adhered to. Conditions of disinfection and rinsing are the responsibility of the physician in charge.

Do not expose the U 8000 S to halogenated solvents (perchloroethylene, chloroform, etc.), halogen containing liquids (e.g. Iodine tinctures), methanol, acetone, phenols or their derivatives.

WARRANTY AND LIMITATION OF LIABILITY

a) The manufacturer warrants that the U 8000 S has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the U 8000 S becoming apparent before the expiry date.

Expiry date: Refer to information on the unit container label.

b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability or fitness or other warranties, which extend beyond those described in paragraph a) above.

The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the U 8000 S and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the U 8000 S, whether as a result of any defect therein or otherwise.

c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings, directions and instructions in the labelling, damage arising from events after the manufacturer's release of the U 8000 S, failure or omission to inspect the U 8000 S before use in order to insure that the U 8000 S is in proper condition, or any warranty given by independent distributors or dealers.

EXPLANATION OF SYMBOLS

= Read instructions before use!

= Expiry date

= Sterilizing agent: Steam

= Lot number

= Recommended storage conditions :
below + 30°C (+ 86F)

= Catalogue no.

= Manufacturer

Directions For Use

CAUTION ! Strict adherence to the procedural order and content, including all operational WARNINGS and CAUTIONS, described in the following procedures by the manufacturer is required to ensure a safe and effective treatment.

CAUTION ! Aseptic technique is required throughout preparation and treatment to avoid contamination.

PRE-TREATMENT

SET UP

The U 8000 S is delivered sterile, with non-pyrogenic fluid compartments. Make certain that the package and the device have not been damaged or exposed to liquid, solvents or other hazardous material.

CAUTION ! Always, check the ultrafilter carefully and use only if it is mechanically intact and the protection caps over the fluid and ultrafiltrate ports are dry and securely in place. If the ultrafilter or packaging is damaged, sterility may be compromised. Therefore, do not use the ultrafilter.

When the U 8000 S is inserted into the incoming water line of a dialysis machine, the pre-ultrafilter pressure must be sufficient to feed the machine used despite the pressure drop in the U 8000 S.

When inserted in the dialysis fluid circuit the U 8000 S is designed for use with a dialysis machine, which controls the ultrafiltration volume/rate. The dialysis machine must be able to compensate for the increased pressure drop in the dialysis fluid circuit caused by the U 8000 S, without interfering with the ultrafiltration control during the treatment session.

This operator's manual describes the use of the U 8000 S when inserted in the dialysis fluid circuit. For operating instructions in the water circuit or concerning the preparation of sterile substitution fluid please contact the dialysis machine manufacturer.

1. Connect the U 8000 S to the dialysis fluid circuit as illustrated in fig. 1. Be certain to use aseptic technique when connecting the sterile filtrate line.

START UP, DEAERATION AND RINSING

2. The U 8000 S shall be rinsed and deaerated during the rinsing of the dialysis fluid circuit.

Thorough deaeration of the U 8000 S during the rinsing procedure is essential for optimal function of the ultrafilter and the dialysis machine.

3. The deaeration shall be done via the deaeration line connected to the deaeration/fluid outlet on the top of the ultrafilter (fig.1). Make sure that there is a proper distribution of fluid through both the lumen and the filtrate compartment during rinsing.

Caution: The U 8000 S must be upright while being deaerated.

4. Clamp the deaeration line from the U 8000 S when the ultrafilter has been thoroughly deaerated.

Caution: If the deaeration line is not clamped, some dialysis fluid will bypass the dialyzer and flow directly to the drain. The dialysis efficiency may be reduced.

TREATMENT

5. When the U 8000 S has been properly rinsed and deaerated, proceed with treatment set up procedures as defined in the respective operator's manuals for the dialysis machine, dialyzer, blood tubing, and clinic protocols.

Caution: An excessively high pressure drop over the U 8000 S might cause problems with the dialysis machine. Problems with an excessively high pressure drop may also lead to a limitation to use low ultrafiltration rates when high flux membranes are involved.

DISINFECTION

Disinfection is possible with the chemical disinfectants listed below or by heat disinfection procedures.

The heat disinfection procedure can include hot citric acid or sodium carbonate according to the concentration stated below.

Disinfection should be carried out after each treatment. When the U 8000 S is used for on-line preparation of non-pyrogenic substitution fluid with the Gambro AK 100/200 ULTRA or AK 200 ULTRA S, the disinfection routine recommended in Operator's manual of the machine should be followed. Rationale: If bacteria are allowed to proliferate undisturbed on the feed side of an ultrafilter, there may be a breakthrough of bacterial and/or endotoxins due to a heavy load on the filter.

Warning ! If contamination of bacteria have occurred on the filtrate side of the ultrafilter or in the tube from the ultrafilter, there will be no possibility to maintain the microbiological quality of the dialysis fluid.

Caution: If the physician in charge decides to chemically disinfect the U 8000 S together with the dialysis machine, the General Warnings and Cautions found in this manual must be observed.

The following chemical disinfection agents are recommended for the U 8000 S:

Peracetic acid:
recommended concentration is according to the manufacturer's instructions, not to exceed 0.1 %

Sodium carbonate
recommended concentration is according to the manufacturer's instructions, not to exceed 2.0 %

Sodium hypochlorite
recommended concentration is according to the manufacturer's instructions, not to exceed 0.5 %

Caution: Chemical disinfectants may be toxic. The user must therefore take necessary precautions before use. Chemical disinfectants may also be harmful to the material used in the fluid path in the dialysis machines.

Caution: The user must validate the disinfection efficiency from a microbiological standpoint.

Note: If the U 8000 S is used in combination with Gambro AK 100/200 ULTRA machine, use only those disinfectants which are recommended in the operator's manual of the AK 100/200 ULTRA. The disinfection procedures should follow those described in the operator's manual of the AK 100/200 ULTRA.

PREPARATION

Caution: Make sure that the chemical disinfection procedure of the machine has not affected the pressure drop caused by the U 8000 S.

1. Open the clamp on the de-aeration line. Proceed with the chemical disinfection procedure according to procedures and instructions valid for the specific machine.

Caution: Ensure that there is a flow of rinsing fluid through both the lumen and filtrate compartment of the U 8000 S during the disinfection procedure.

2. After the disinfectant dwell time, specified in the operator's manual for the dialysis machine, the system must be rinsed to eliminate the disinfectant.

Caution: The user must validate the disinfection efficiency from a microbiological standpoint.

Caution: Chemical disinfectants are toxic. The residual concentration in the dialysis fluid must therefore be below the levels specified by the clinic or by national standards. It is essential to use a test method recommended by the manufacturer of the disinfectant.

DECALCIFICATION

The ultrafilter U 8000 S can be included in the decalcification procedure with hot citric acid. The citric acid concentration must not exceed 2%.

General Information

PERFORMANCE

The performance values given below should be regarded as approximate. When the ultrafilter is used under clinical conditions, different values may be obtained due to the clinical setting, measuring techniques, and variations in manufacturing.

PERFORMANCE		U 8000 S	
Determined with 0,9% sodium chloride solution. Mean pressure drop is measured between fluid inlet and filtrate outlet.			
Filtrate flow Q_F (ml/min)	Pressure drop. Δp (mmHg)		
	solution 21°C		solution 37°C
300	43 - 64	29 - 43	
500	71 - 107	48 - 71	
700	100 - 150	67 - 100	

Retention of bacteria	
Type of challenge	LRV *
Bacterial: <i>Pseudomonas diminuta</i> ATCC 19146	> 7
Endotoxin	
Typical value for endotoxin after filtration:	< 0,05 IU/ml
* LRV = Logarithmic reduction value;	
$LRV = \log_{10} \frac{\text{number of organisms in challenge suspension}}{\text{number of organisms in filtrate}}$	

Upon request, information about methods used to obtain performance data is available.


Figure 1: Recommended connection of the U 8000 S when used in the dialysis circuit.
(Depending on the type of dialysis machine, the connections may vary)

SPECIFICATIONS

U 8000 S	
Effective membrane area (m ²)	2.1
Fiber dimensions	
Inner diameter (µm)	215
Wall thickness (µm)	50
Priming volume (ml)	
lumen	152
filtrate side	306
Maximum TMP (mmHg)	600
Sterilization agent	Steam
Sterile barrier	Medical grade paper
Components	Materials
Membrane	Polyamix™ (PAES/PVP/PA)
Potting material	Polyurethane (PUR)
Housing and caps	Polycarbonate (PC)
Sterile plugs	Polypropylene (PP)

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