

产品名 GamCath 留置导管

适用范围 GamCath 留置导管适于短期应用。它可用于建立血管通道以进行体外血液净化。它可用于锁骨下静脉、颈静脉或股静脉置管。

型号、规格	产品名称	订购号	外径(Fr)	长度(cm)	其他信息	套/盒
	GDC-1115	101797	11	15	-	5
	GDC-1117.5	101242	11	17.5	-	5
	GDC-1120	101799	11	20	-	5
	GDC-1125	101800	11	25	-	5
	GDC-1112.5J	101028	11	12.5	弯头	5
	GDC-1115J	101798	11	15	弯头	5
	GDC-1120J	101295	11	20	弯头	5
	GDK-1115	955055	11	15	-	5
	GDK-1117.5	103241	11	17.5	-	5
	GDK-1120	955056	11	20	-	5
	GDK-1125	101511	11	25	-	5
	GDK-1112.5J	955053	11	12.5	弯头	5
	GDK-1115J	955054	11	15	弯头	5
	GDK-1117.5J	101379	11	17.5	弯头	5
	GDK-1120J	955057	11	20	弯头	5
	GDK-1215	100842	12	15	-	5
	GDK-1220	101803	12	20	-	5
	GDK-610P	100943	6	10	-	5
	GDK-612.5P	955059	6	12.5	-	5
	GDK-810P	101512	8	10	-	5
	GDK-812.5P	955058	8	12.5	-	5
	GTK-1215	101808	12	15	-	5
	GTK-1217.5	103482	12	17.5	-	5
	GTK-1220	101809	12	20	-	5
	GTK-1215J	101384	12	15	弯头	5
	GTK-1220J	101385	12	20	弯头	5
	GDHK-1115	955424	11	15	-	5
	GDHK-1120	955425	11	20	-	5
	GDHK-1125	106661	11	25	-	5
	GDHK-1315	955060	13	15	-	5
	GDHK-1317.5	103357	13	17.5	-	5
	GDHK-1320	955061	13	20	-	5
	GDHK-1325	955062	13	25	-	5
	GDHK-1115J	112572	11	15	弯头	5

GDHK-1120J	112573	11	20	弯头	5
GDHK-1315J	109535	13	15	弯头	5
GDHK-1320J	109511	13	20	弯头	5
MS-GDHK-1315	955063	13	15	-	5
MS-GDHK-1317.5	103521	13	17.5	-	5
MS-GDHK-1320	955064	13	20	-	5
MS-GDHK-1325	106662	13	25	-	5

产品图

